

## ASP-18

Components and Final Concentration in Culture Medium	Stock Solution	Addition per Litre of Culture Medium
1. HEPES (3.00 mM)	238.10 g / l dH <sub>2</sub> O	3 ml
2. NaCl (479.00 mM)		weigh and add 18 g
3. KCl (9.40 mM)	60.00 g / l dH <sub>2</sub> O	11.6 ml
4. MgSO <sub>4</sub> x 7 H <sub>2</sub> O (28.40 mM)		weigh and add 7 g
5. MgCl <sub>2</sub> x 6 H <sub>2</sub> O (19.70 mM)		weigh and add 4 g
6. CaCl <sub>2</sub> x 2 H <sub>2</sub> O (10.00 mM)	370.00 g / l dH <sub>2</sub> O	4 ml
7. NaNO <sub>3</sub> (1.18 mM)	100.30 g / l dH <sub>2</sub> O	1 ml
8. K <sub>3</sub> PO <sub>4</sub> x 3 H <sub>2</sub> O (47.00 μM)	12.50 g / l dH <sub>2</sub> O	1 ml
9. Na <sub>2</sub> -Glycerophosphate (31.70 μM)	6.85 g / l dH <sub>2</sub> O	1 ml
10. Na <sub>2</sub> SiO <sub>3</sub> x 9 H <sub>2</sub> O (528.00 μM)	28.42 g / l dH <sub>2</sub> O	5.3 ml
11. NTA (Titriplex I) (523.00 μM)	10.00 g / l dH <sub>2</sub> O	10 ml
<b>12. Vitamin Solution</b>		1 ml
Vitamin B12 (0.15 nM)	0.20 mg / l dH <sub>2</sub> O	
Biotin (4.10 nM)	1.00 mg / l dH <sub>2</sub> O	
Thiamine-HCl (0.30 μM)	100.00 mg / l dH <sub>2</sub> O	
Niacinamide (0.80 nM)	0.10 mg / l dH <sub>2</sub> O	

pH of the Vitamin Solution should be around pH 7.0

### 13. Trace Metals

Na<sub>2</sub>EDTA x 2 H<sub>2</sub>O: 4.36 g

FeCl<sub>3</sub> x 6 H<sub>2</sub>O: 3.15 g

Dissolve in 1000 ml dH<sub>2</sub>O, then add 1 ml of Primary Trace Metals each (see below).

Primary Trace Metals are stored frozen as 1 ml aliquots.

#### 13.2. Primary Trace Metals

13.2.1. K <sub>2</sub> CrO <sub>4</sub>	0.194 g / 100 ml dH <sub>2</sub> O
13.2.2. CoCl <sub>2</sub> x 6 H <sub>2</sub> O	1.00 g / 100 ml dH <sub>2</sub> O
13.2.3. CuSO <sub>4</sub> x 5 H <sub>2</sub> O	0.25 g / 100 ml dH <sub>2</sub> O
13.2.4. MnCl <sub>2</sub> x 4 H <sub>2</sub> O	18.00 g / 100 ml dH <sub>2</sub> O
13.2.5. Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	1.89 g / 100 ml dH <sub>2</sub> O
13.2.6. NiSO <sub>4</sub> x 6 H <sub>2</sub> O	0.27 g / 100 ml dH <sub>2</sub> O
13.2.7. H <sub>2</sub> SeO <sub>3</sub>	0.13 g / 100 ml dH <sub>2</sub> O
13.2.8. Na <sub>3</sub> VO <sub>4</sub>	0.184 g / 100 ml dH <sub>2</sub> O
13.2.9. ZnSO <sub>4</sub> x 7 H <sub>2</sub> O	2.20 g / 100 ml dH <sub>2</sub> O

### 14. Soil Extract (optional)

See below

10 ml

adjust the pH to 7,9 and autoclave

#### Preparation of Soil Extract

10 g of garden-soil is mixed with 120 ml dH<sub>2</sub>O and boiled for 10 minutes. Afterwards it is centrifuged for 10 minutes (low speed), and the supernatant is filtered through a series of membrane filters from 1.2 μm – 0.1 μm pore size. The remaining filtrate is adjusted to 100 ml with dH<sub>2</sub>O. Aliquots of 10 ml are stored frozen.

The soil should not be recently fertilized and should not contain too much humus.