

## SFM + Si

Components and Final Concentration in Culture Medium	Stock Solution	Addition per Litre of Culture Medium
1. HEPES (1.00 mM)	238.10 g / l dH <sub>2</sub> O	1.0 ml
2. Ca(NO <sub>3</sub> ) <sub>2</sub> x 4 H <sub>2</sub> O (0.21 mM)	100.00 g / l dH <sub>2</sub> O	0.5 ml
3. MgSO <sub>4</sub> x 7 H <sub>2</sub> O (0.203 mM)	20.00 g / l dH <sub>2</sub> O	2.5 ml
4. K <sub>2</sub> HPO <sub>4</sub> x 3 H <sub>2</sub> O (13.20 μM)	5.00 g / l dH <sub>2</sub> O	0.6 ml
NaNO <sub>3</sub> (0.35 mM)	50.00 g / l dH <sub>2</sub> O	
Na <sub>2</sub> CO <sub>3</sub> (0.19 mM)	32.00 g / l dH <sub>2</sub> O	
5. H <sub>3</sub> BO <sub>3</sub> (16 μm)	1.00 g / l dH <sub>2</sub> O	1 ml
<b>6. Vitamin Solution</b>		1 ml
Vitamin B12 (0.15 nM)	0.20 mg / l dH <sub>2</sub> O	
Biotin (4.10 nM)	1.00 mg / l dH <sub>2</sub> O	
Thiamine-HCl (0.30 μM)	100.00 mg / l dH <sub>2</sub> O	
Niacinamide (0.80 nM)	0.10 mg / l dH <sub>2</sub> O	
pH of the Vitamin Solution should be around pH 7		
<b>7. Trace Metals</b>		1 ml
7.1. Preparation of Trace Metal Solution		
Na <sub>2</sub> EDTA x 2 H <sub>2</sub> O : 4.36 g		
FeCl <sub>3</sub> x 6 H <sub>2</sub> O : 3.15 g		
Dissolve in 1000 ml dH <sub>2</sub> O, then add 1 ml of Primary Trace Metals each (see below).		
Primary Trace Metals are stored frozen as 1 ml aliquots.		
<b>7.2. Primary Trace Metals</b>		
7.2.1. K <sub>2</sub> CrO <sub>4</sub>	0.194 g / 100 ml dH <sub>2</sub> O	
7.2.2. CoCl <sub>2</sub> x 6 H <sub>2</sub> O	1.00 g / 100 ml dH <sub>2</sub> O	
7.2.3. CuSO <sub>4</sub> x 5 H <sub>2</sub> O	0.25 g / 100 ml dH <sub>2</sub> O	
7.2.4. MnCl <sub>2</sub> x 4 H <sub>2</sub> O	18.00 g / 100 ml dH <sub>2</sub> O	
7.2.5. Na <sub>2</sub> MoO <sub>4</sub> x 2 H <sub>2</sub> O	1.89 g / 100 ml dH <sub>2</sub> O	
7.2.6. NiSO <sub>4</sub> x 6 H <sub>2</sub> O	0.27 g / 100 ml dH <sub>2</sub> O	
7.2.7. H <sub>2</sub> SeO <sub>3</sub>	0.13 g / 100 ml dH <sub>2</sub> O	
7.2.8. Na <sub>3</sub> VO <sub>4</sub>	0.184 g / 100 ml dH <sub>2</sub> O	
7.2.9. ZnSO <sub>4</sub> x 7 H <sub>2</sub> O	2.20 g / 100 ml dH <sub>2</sub> O	
8. Na <sub>2</sub> SiO <sub>3</sub> x 9 H <sub>2</sub> O (0.50 mM)	28.42 g / l dH <sub>2</sub> O	5 ml

adjust the pH to 7,0 and autoclave