


Standard Operating Procedure
**Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

Starting Up the System

- Start up the computer, start Windows and log in as Operator
- Turn on the cytometer main power
- Start BD FACSDiva software with your Group-name
- Turn on the Laser power (for best performance wait 30' for laser stabilization)
- Check fluid levels in the Cytometer window





Fluidics Startup

- Go to the task bar and choose "cytometer", "fluidics startup"
(If there is air in the sheath filter: take a syringe, plug it into the wheel, screw it open and suck up liquids and air until there is no air anymore, close the screw and)
- Follow the program guide (but take off the closed loop nozzle from the flow cell)

Special setup

- Switch on both switches of the water bath-temperature unit
 - Press the power button of the temperature unit
 - Make sure there is enough distilled water in the chamber
 - Make sure it's set up for your desired temperature
 - Connect input and output tubing to the tube holder
 - Press the power button



- For **S2 sorting** switch on the air pressure system under the table.
 - Check the life-time (if only red light is left - tell the staff)
 - Make sure there is not more than 20% of pressure applied



**Standard Operating Procedure
Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

- Make sure the deflection plates are clean (important for sorting!)
- Verify that the appropriate sort setup is selected. (Go to Cytometer -> View Configuration ; choose the appropriate nozzle setup in the new window; press ->set configuration->OK and close the window)
- Turn on the stream
- Optimize the breakoff
 - with the Amplitude (ONLY)
 - Break off point needs to be in the upper third
 - Max. 5 satellite drops
 - Watch out for the correct gap-size (nozzle size dependent)

Nozzle	Cell size	Gap
70µm	<12µm	6-7
85µm	>12µm	8-10
100µm	≤15µm	10-12



Now transfer "drop 1" and "gap" actual-value in the stream window to the white target-value field left beside it.

- Turn **on** the Sweet Spot
- Prepare a tube with one Drop of CS&T beads with ~400µl PBS
- Use Experiment ->New experiment -> CS&T test tube and run as a sample
 - OR Perform CST system setup in Cytometer ->CST(this will take 15 Minutes)
 - Take care to setup the correct CST lot number
 - Run the performance check

Ampl:	24.4
Freq:	87.0
Drop 1:	198 198
Gap:	6 8

Target value Actual value

**Standard Operating Procedure
Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

Adjust drop delay (every experiment)


- Use 1000µl PBS with 1 drop of “Accu Drop Beads” from the fridge
- Go to task bar and use “new experiment” choose “Accu Drop delay”
- Klick on “+” before the “specimen”, “tube”,
- Activate the little arrow on the left side (tube pointer) of the “tube_001” so it turns green
- doubleclick on “sort layout”
- **Close the cover of the machine**
- Load the tube and acquire
- Adjust the “threshold rate” with the “flow rate” in the acquisition dashboard.
- If you need to increase the “flow rate” higher than 5 please dilute the accudrops!

Nozzle size	Flow rate events
70nm	1000-3000
85nm	1000-2000
100nm	600-1500

- Switch on “voltage” and the “test sort” –button
- Activate the **near-left** side-stream
- Use the silver screw at the sorting chamber in the machine to get the stream-dots (most important in this case is the middle stream dot and the near left one.) as bright as possible
- Activate “optical filter” button;
 - to make sure to fit the near-left side stream in the filter box just deactivate and activate the optical filter
- **Close the cover of the machine**
- start the sort (**cancel !** warning window)
- Click on “auto delay”
- Use “start run”-button
 - => ! the result should be a parabolic curve with its maximum in the middle !
- Unload the tube
- Don´t save a copy of the sort report

Adjust the stream

- Put two or four tubes in the tube holder; adjust the sidestreams with the voltage sliders roughly
- Open machine, open sorting chamber;
- Turn “**voltage**” on again
- **BE CAREFUL NOW! DON´T TOUCH THE DEFLECTION PLATES!**
- Choose “test sort”, close “waste drawer” and
- If the streams don´t hit the tubes properly adjust them with the voltage sliders.
- Turn off “test sort”, “voltage”, “waste drawer”, close the chamber and the machine

Imaging Center Essen	Institut für Experimentelle Immunologie und Bildgebung	 Universitätsklinikum Essen
<p style="text-align: center;">Standard Operating Procedure Workflow BD FACS Aria III Cell Sorter</p> <p><u>Please note:</u></p> <p>If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.</p>		<p style="text-align: right;">Page 4 von 10</p> <p style="text-align: right;">Last change: 25.02.2020</p>


Setting Up an experiment with compensation

- Create a new experiment and use e.g. the Doublet Discrimination analysis template.
- Select “cytometer settings” and tick height and width for the FSC and SSC parameters in the cytometer window.

- Delete all colors you do not need!


- Select “Experiment” > “Compensation Setup” > “Create Compensation Controls” > “OK”
- Select “compensation control”; double-click on “unstained”
- Load the unstained control tube.
- Set voltages for FSC and SSC appropriate to your cell type
- Adjust a gate around the population of interest, right click and “apply to all compensation controls”
- Load all tubes in the unstained control to setup the voltages properly.
- Start again with the unstained sample and “record data”
- Now record the data for each tube correctly.
- Check your single staining and adjust the bar around your brightest positive population if not done by machine properly after recording.
- Finish the compensation
- “Select Experiment” > “Compensation Setup” > “Calculate Compensation”
- If everything was correct choose “apply” or “**link and save**” to save your compensation to the panel
- Check your settings in “cytometer” “compensation “ if values are taken over

- You might want to perform (Fluorescence minus one) “ FMO-controls” to confirm that your compensation is correct

Imaging Center Essen	Institut für Experimentelle Immunologie und Bildgebung	 Universitätsklinikum Essen
<p style="text-align: center;">Standard Operating Procedure Workflow BD FACS Aria III Cell Sorter</p> <p><u>Please note:</u></p> <p>If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.</p>		<p style="text-align: right;">Page 5 von 10</p> <p style="text-align: right;">Last change: 25.02.2020</p>

Sorting without additional compensation

- Create a new experiment and apply e.g. the Doublet Discrimination analysis template or use an existing experiment.
- Select “cytometer settings” and tick height and width for the FSC and SSC parameters in the cytometer window.
- Delete colors you do not need
- Create a “new tube” and apply a “sort layout” if not there yet
- Populate the global worksheet with plots (dot plot / density plot / histogram)
- Set all your gates of interest with your unstained cells
- Setting up your gating strategy:
 - gates for all tubes in your experiment on a global worksheet
 - gates for specific, single tubes in tube specific worksheet
- acquire some data
 - You might want to perform (Fluorescence minus one) FMO-controls to confirm that your manual settings are correct
- To create further tubes with the same settings just use the “net tube” button beside the “load” button
- Choose which population should go to which tube in the sort layout
 - ☞ big cells and populations in the middle, the smaller ones in the outer tubes
- If you apply stopping criteria do this for your smallest population
- Install the collection tubes and check the side streams
- If necessary change collection tubes and **fill in some sample buffer** (1 sorting droplet is ~1nl)
- “Load” your sample, “record”, “sort”
- **Don't cancel** the warning window
- Monitor the stream during sorting.
- Stop the sort and save a copy of the sort report if needed
- Take away the tube holder when you are finished

Imaging Center Essen	Institut für Experimentelle Immunologie und Bildgebung	 Universitätsklinikum Essen
<p style="text-align: center;">Standard Operating Procedure Workflow BD FACS Aria III Cell Sorter</p> <p><u>Please note:</u></p> <p>If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.</p>		<p style="text-align: right;">Page 6 von 10</p> <p style="text-align: right;">Last change: 25.02.2020</p>


➤ **If there is a user after you :**
Standby

- Run filtered **FACS rinse** ~5 minutes FR11 (Flow Rate 11)
- Run filtered **FACS clean** for ~5 minutes FR11
- Run filtered **dH₂O** ~5 minutes FR11
- Logout (yes...let the stream keep on running)
- Open the cover of the machine

➤ **If you are the last user:**
Shutting down the system

- Run filtered **FACS rinse** for ~2 minutes FR11 (Flow Rate 11)
- Run filtered **FACS clean** ~2 minutes FR11
- Run filtered **dH₂O** ~2 minutes FR11
- Turn **off** the stream
- Change to the closed loop nozzle
 - Then put a full tube of FACS Rinse in to the loading port
 - Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell”
 - Put a full tube of FACS Clean in to the loading port
 - Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell”
 - Put a full tube of dH₂O in to the loading port
 - Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell” twice
- Go to “cytometer” in the task bar and **perform a “fluidics shutdown” on Friday**
 - Follow the wizard
 - The wizard asks for “Cleaning fluidics”. This is **dH₂O**
- Turn off the Laser power
- Turn off the cytometer main power
- Shut down the computer
- Shut down the cooling unit
- Switch back the connectors from EtOH to PBS tank

Fill out the users-template in the folder.

Imaging Center Essen	Institut für Experimentelle Immunologie und Bildgebung	 Universitätsklinikum Essen
<p style="text-align: center;">Standard Operating Procedure Workflow BD FACS Aria III Cell Sorter</p> <p><u>Please note:</u></p> <p>If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.</p>		<p style="text-align: right;">Page 7 von 10</p> <p style="text-align: right;">Last change: 25.02.2020</p>

If you want to sort S2 specimen, or human blood samples please follow the IMCES S2 working procedure!

Before sorting S2-cells please:

- Switch on the air pressure system under the table.
 - Check the filter- life-time (if only red light is left - tell the staff)
 - Make sure there is not more than 20% of pressure applied
- For your own safety make sure you touch your sample, as well as the Aria only with gloves.

After sorting either human , or S2-cells please:

- Run filtered **FACS rinse** for ~2 minutes FR11 (Flow Rate 11)
- Run filtered **FACS clean** ~2 minutes FR11
- Run filtered **dH₂O** ~2 minutes FR11

- Stop the stream
- Change to the closed loop nozzle
- Then put a **full** tube of FACS Rinse in to the loading port
- Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell”
- Put a **full** tube of FACS **Clean at 40 °C!!!** in to the loading port
- Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell”
- **Leave this hot solution in the flow cell for 2 minutes**
- Put a **full** tube of ddWater in to the loading port
- Go to the task bar and choose: “cytometer”, “cleaning modes”, “clean flow cell” twice

- Go to “cytometer” in the task bar and **perform a “fluidics shutdown”**
 - Follow the wizard
 - The wizard asks for “Cleaning fluidics”. This is **dH₂O**

- Turn off the Laser powers
- Turn off the cytometer main power
- Shut down the computer
- Turn off the cooling unit
- Switch back the connectors from EtOH to PBS tank

**After your sort clean the Aria and your bench space and the computer table
with 70% Ethanol!**



**Standard Operating Procedure
Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

Troubleshooting

Problem:	Precise Problem:	Possible solutions:
Stream; there is no stream	<ul style="list-style-type: none"> There is no visible stream in the stream chamber 	<ul style="list-style-type: none"> Activate the stream in the stream window Check if there is a nozzle in the flow cell Check if the pressure tubing is connected to the PBS tank Check if the blue PBS tube is connected to the blue tubed filter Check if there is enough PBS in the tank
	<ul style="list-style-type: none"> There is a stream in the stream chamber but not visible in the Software 	<ul style="list-style-type: none"> Clean the camera besides the stream
Stream; stream is not looking as usual	<ul style="list-style-type: none"> Stream is vibrating 	<ul style="list-style-type: none"> Check if there is the tank vibrating against the blue tubed filter Check if something is hitting along the blue tube
	<ul style="list-style-type: none"> Stream is not forming droplets 	<ul style="list-style-type: none"> Make sure the right sort setup is chosen for the inserted nozzle Check if the right nozzle is inserted Make sure the blue tube is not connected to the ethanol tank
Accu Drops; there are no stream dots	<ul style="list-style-type: none"> There are no visible side streams in the stream chamber 	<ul style="list-style-type: none"> Make sure you are sorting with PBS Check if the voltage is switched on Activate the test sort button again Make sure the side stream voltage sliders are >0
	<ul style="list-style-type: none"> There is just a bright line in the side stream window 	<ul style="list-style-type: none"> The stream might be hitting the waste; adjust the chamber with the screw driver



**Standard Operating Procedure
Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

	<ul style="list-style-type: none"> • There are side streams visible in the stream chamber but not in the software 	<ul style="list-style-type: none"> • Close the cover of the machine if its open • Use the silver screw to bring the stream laser in the focal plane
Accu Drops; the side streams don't look as usual	<ul style="list-style-type: none"> • The center stream does not look as usual 	<ul style="list-style-type: none"> • Use the silver screw to bring the stream laser in the focal plane • Use the 2nd 3rd 4th Drop values to tighten the center stream dot (every value should be half of the one before) • Change the Amplitude by +/- 1
	<ul style="list-style-type: none"> • You cannot get the center stream in the "optical filter" 	<ul style="list-style-type: none"> • Use the center stream voltage slider to adjust • Move the camera behind the stream chamber with your hand
	<ul style="list-style-type: none"> • The side streams don't look as usual 	<ul style="list-style-type: none"> • Use the silver screw to bring the stream laser in the focal plane
	<ul style="list-style-type: none"> • There are more side streams than desired 	<ul style="list-style-type: none"> • Make sure there are only the side streams activated you wanted (set all other sider stream voltage sliders to 0) • Change the 2nd 3rd 4th Drop values (every value should be half of the one before)
	<ul style="list-style-type: none"> • You cannot see the left stream 	<ul style="list-style-type: none"> • Check if the optcal filter is activated; deactivate • Make sure the side stream voltage slider is not set to 0



**Standard Operating Procedure
Workflow BD FACS Aria III
Cell Sorter**

Please note:

If images acquired on IMCES instruments are used in publications we would be grateful if you could mention the facility in the Acknowledgements.

<p>Accu Drops; the Accu Drop Delay is <u>not</u> working</p>	<ul style="list-style-type: none"> The machine is not measuring beads after loading the tube 	<ul style="list-style-type: none"> Make sure the cover of the machine is closed Open an new Accu Drop Delay- template via the task bar Check if the lasers are switched on
	<ul style="list-style-type: none"> You cannot start the measurement 	<ul style="list-style-type: none"> Make sure the sorting is not paused Check if you started the Auto Delay Adjust the maximum threshold rate for your nozzle
	<ul style="list-style-type: none"> The Accu Drop Delay is not coming to a result 	<ul style="list-style-type: none"> Check if you set up the near left side stream and not the far left
	<ul style="list-style-type: none"> There is no curve appearing after measurement 	<ul style="list-style-type: none"> Just compare the result to a prior measurement with the same nozzle and the same Frequency
	<ul style="list-style-type: none"> The curve is not parabolic 	<ul style="list-style-type: none"> Check if you are measuring Accu Drop Beads Make sure the deflection plates are clean Rerun the Accu Drop Delay