

How to dispose syringes and needles

About the right way to handle them

Needles

- It's not allowed to refit the cap to the used needle
- injury risk!



Needles

- Needles must not be disposed into the laboratory waste
- injury risk!



needles

- Used needles have to be disposed directly at the workplace into the dedicated containers
- The container has a specific appliance for pulling off the needle from the injection syringes



Needles

- These containers are only for needles!
- Not for paper, plastic caps or syringes!



Syringes

- Disposable syringes have to be disposed into the lab waste
- Except biohazard contaminated syringes, which are to collect into designated containers to autoclave or tyndallize them!