

Euglena media

Components are dissolved in 950 ml dH₂O and filled up to 1 liter.

Natriumacetat x 3 H₂O		1 g
Beef Extract		1 g
Pepton		2 g
Yeast Extract		2 g
Soilextract		10 ml
Vitamin Solution		1 ml
Vitamin B12 (0.15 nM)	0.20 mg / l dH ₂ O	
Biotin (4.10 nM)	1.00 mg / l dH ₂ O	
Thiamine-HCl (0.30 μM)	100.00 mg / l dH ₂ O	
Niacinamide (0.80 nM)	0.10 mg / l dH ₂ O	

pH of the Vitamin Solution should be around pH 7